PURITY

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PURITY

Operating Manual

It is essential that these notes should be carefully read before operating the unit since an incorrect use could affect the correct operation of the sterilizer.

In case of doubts do not hesitate to contact the reseller that supplied the product.

Thank you for your trust.

All rights reserved.

No part of this publication can be reproduced, transmitted, transcribed, filed in an information retrieval system or translated in any language, or computer languages, in any form or media, without previous written authorization by the company **Dental X**.

Dental X reserves the right to make changes in the specifications of the product described in this manual at any time without previous notice and without the inform anyone about these changes or revisions.

PURITY is a registered trademark of Dental X.

STERILINE is a registered trademark of Dental X.

DENTAL X S.p.A.

Via Marzotto 11

36031 Dueville (VI) Italy

tel. +39 0444 367400

Fax +39 0444 367436

e-mail: dentalx@dentalx.it

<http://www.dentalx.it>

TECHNICAL SUPPORT

Italy - tel. +39 0444 367415

Outside Italy - tel. +39 02 76275234

CONTENTS

GENERAL INFORMATION	1
INTRODUCTION.....	1
CONFORMITY TO EUROPEAN DIRECTIVES	1
1. FIRST LOOK	2
1.1 PACKAGE SIZE AND WEIGHT	2
1.2 EXTERNAL DIMENSIONS AND USEFUL MEASUREMENTS.....	2
1.3 SAFETY FEATURES.....	2
1.4 PACKAGING CONTENT	3
1.5 TECHNICAL FEATURES	3
1.5.1 Environmental operating conditions.....	3
1.6 INDICATION OF USE	3
2. INSTALLATION	4
2.1 BASIC REQUIREMENTS	4
2.2 ASSEMBLY.....	4
2.3 CONNECTIONS	4
2.4 FIRST START-UP	5
3. USE INSTRUCTIONS	6
3.1 COMMANDS/SIGNALS.....	6
3.2 OPERATING.....	6
4. MAINTENANCE	7
4.1 GENERALS	7
4.2 CLEANING OF THE INPUT PRE-FILTER	7
4.3 DIAGNOSTICS.....	7
4.4 FILTER KIT REPLACEMENT.....	8

GENERAL INFORMATION

INTRODUCTION

This manual provides instructions for:

- correct installation
- correct use
- full maintenance

The unit should be installed following the instructions in this manual.

Please note that the customer is responsible for the fulfillment of the legal concerning the installation and operation of this unit.

Should the unit not be properly installed and operated, as well as in case of inadequate maintenance, the manufacturer cannot be liable for any damage, injury or incorrect operation.

We recommend that you first check that the package is intact and that the included components match the contents of the delivery bill.

NOTE **In case of apparent damage or missing items, immediately inform the carrier and send notice to the Dental X company or its area representative.**

CONFORMITY TO EUROPEAN DIRECTIVES

The water processing system PURITY manufactured by DENTAL X conforms to the EN 1717 norm as concerns the protection against reflux towards the main waterline.

This product by Dental X has been designed and manufactured using high quality materials and components that can be recycled and reused.




This symbol means that electrical and electronic appliances, at the end of their life cycle, must be disposed separately from domestic waste; dispose this unit by delivering it to a collection/recycling center. Legal sanctions are applicable if the regulation concerning waste disposal should not be met.

The used filter cartridges can be disposed as simple domestic waste.


Help us to preserve the environment in which we live!

1.FIRST LOOK

1.1 PACKAGE SIZE AND WEIGHT

Rear view

Overall weight: 10 Kg when dry, 12 Kg in use
Keep the original packing.

1.2 EXTERNAL DIMENSIONS AND USEFUL MEASUREMENTS

Length = 450mm

Depth = 510mm

Overall height = 220mm – Height of the installed unit = 140mm

1.3 SAFETY FEATURES

The unit operates at low voltage and is equipped with an dual insulation external power supply module, protected against over temperatures and voltages.

To ensure safety continuity also in case of fault, replace the power supply module with an original spare part.

PURITY has been designed to avoid possible leaks by collecting the water into the sterilizer's internal tank, where the water in excess is automatically discharged directly into the drain.

Even though the unit is equipped with these protections, it should be switched off when not in use and the water supply tap should be closed when the unit is not attended.

1.4 PACKAGING CONTENT

- PURITY unit
- Power supply module
- PURITY feeding pipe, 1.5 m long, 3/4" connector
- PURITY draining pipe, 2 m long, 6/4 mm across
- Sterilizer draining pipe, 2 m long , 10/14 mm across
- Operating manual
- Warranty card

NOTE For the warranty to be valid, it is required that a copy of the warranty sheet should be sent, through the reseller, to the manufacturer, on pain of warranty voidance.

1.5 TECHNICAL FEATURES

Power supply:	230Vac 50Hz / 24Vac - 800mA
	120Vac 60Hz / 24Vac - 800mA
Max electric consumption:	25W
Water supply:	0,4-6 bar, min. 2 l/minute
Max water consumption:	2 l/minute
Drain:	max 2 l/minute
Water efficiency: 20-40%:	the water draining from PURITY, which is separately directed from the sterilizer draining water, can be recycled for non-alimentary purposes.

1.5.1 Environmental operating conditions

The unit has been designed to operate in environment with temperature range between 3°C and 40°C, relative humidity lower than 95%, pressure from 750 mBar to 1050 mBar and altitude between 0 and +2500 meters above water level

DO NOT USE THE EQUIPMENT IN PRESENCE OF ANESTHETIC FLAMMABLE GAS.

1.6 INDICATION OF USE

PURITY has been designed to feed with demineralized water, obtained through inverted osmosis, the sterilizers manufactured by DENTAL X and properly configured for this feature. The sterilizer compatibility with the PURITY unit is pointed out by a square-shaped plastic cap on the main tank.

2. INSTALLATION

2.1 BASIC REQUIREMENTS

Check that the voltage of your installation matches the value indicated by the label on the power supply module, and that the socket to which the unit will be plugged has a rating of at least 1A.


The manufacturer is not responsible for damage to people or things caused by inadequate electrical installations.

2.2 ASSEMBLY

Remove the sterilizer's tank cover and store it for possible future use.

Lay PURITY on the sterilizer, checking that the three guides fit in the tank without damaging it and that the float is free to move.


WARNING The float is fragile, handle with care.

2.3 CONNECTIONS

Connect the supplied PURITY draining pipe (the pipe with a shorter section) to the central connector on the PURITY back panel; tighten the fixing ring nut. Connect the other end of the pipe to a drain outlet.

Note: *The quantity of drain water is minimal, therefore the pipe can also be let to drain freely in a basin or be directed to the house water drain.*

WARNING The draining pipe is under pressure, and possible noise should be generated if connected directly into a drain; install a discharge aerator.

Local regulation could require the installation of drain separating element.

Use the supplied feeding pipe; connect the elbow fitting end to the connector on the left of the PURITY back panel and the other end to the tap water outlet.

WARNING Provide a safety water tap on the outlet of the water supply.

Use the supplied sterilizer draining pipe; connect the quick coupling fitting to the connector on the sterilizer rear panel (see operating manual for more details) and connect the other end to a draining outlet.

The water draining from the sterilizer only occurs by the effect of gravity; for a correct operation therefore, a minimum inclination of the sterilizer should be arranged together with enough inclination of the draining pipe to allow the water to flow.

Draining water can be hot!.

Connect the output cable plug of the power supply module to the jack socket on the PURITY back panel and the mains plug to an electrical socket.

2.4 FIRST START-UP

WARNING These operations should be carried out only by qualified technicians. Incorrect operations could affect the quality of the sterilizing process.

1. Open the water feeding tap and check that there are no leaks on the connections.
2. Drain the sterilizer's clean water tank.
3. Switch on PURITY by means of the front panel switch (position I).
After a short time, the water will start to flow from the PURITY draining pipe and, after 1-3 minutes, from the nozzle placed above the tank of the sterilizer.
To fill up the sterilizer tank requires 10-20 minutes.

WARNING Discard the first water produced; its quality might not be adequate. Switch off PURITY and drain the sterilizer's clean water tank again.

4. Switch PURITY back on.

WARNING Do not leave the unit on for long periods without a supply of water; the pump could be damaged.

5. Check that the sterilizer draining pipe is free and without obstructions; the turning on of Led MAX means that the maximum level of water in the sterilizer's recycling tank has been reached, and consequently that the draining connection is not working properly.
6. Take note of the installation date and of the number of cycles performed by the sterilizer.

NOTE The filters of the PURITY equipment cannot be tested with water in the internal circuits during the assembling process, as this could cause a decay in performance; the check of the hydraulic tightness is therefore carried out with compressed air, and it is possible that this procedure does not reveal small leaks. We recommend that, at the end of the installation, the PURITY upper cover (fixed with two screws on each side) be removed in order to check for internal water leaks.

3. USE INSTRUCTIONS

3.1 COMMANDS/SIGNALS

The PURITY front panel is provided with the power supply switch; the switch can be left to ON for all the period during which the sterilizer is operative. Electrical consumption is limited to the feeding phase. At the end of the day, switch PURITY off and close the water supply tap.

The signaling of the max water level in the clean water tank appears on the sterilizer front panel.

The sterilizer will signal, through the internal check of the water quality, when the PURITY filters are used up and should therefore be replaced.

3.2 OPERATING

The PURITY unit is composed by serial connected filters: mechanical pre-filter, pump, activated carbon filter, reverse osmosis membrane and ion-exchange resin filter.

The membrane performs the separation of polluting chemicals that are sent to the draining pipe; therefore the draining water has a concentration in polluting chemicals higher than the input water.

The float in the sterilizer's clean water tank controls the water level and drives PURITY when needed.

The sterilizer's water tank has a capacity for several sterilization cycles (the number depends on the model and on the type of cycle), therefore it is not necessary for PURITY to be always on; however, keep in mind that, to make available the minimum quantity of water required for a cycle, several minutes are required.

The water used during a sterilization cycle is condensed and normally collected in the sterilizer's recovery tank. When the water level in the tank exceeds half its capacity, the used water spills over in the draining pipe; in this way it will no longer necessary to fill the main tank nor empty the recovery one.

Anyway, if needed, the sterilizer can be normally used without connecting the PURITY unit.

4. MAINTENANCE

4.1 GENERALS

The PURITY unit does not use mechanical filtering alone, as in a simple demineralizer system. The PURITY lifetime is therefore much longer and the unit does not require special maintenance until it has provided 1000-2000 liters of water.

The maintenance period depends on the quality of the supply water; the experience on the field indicates that, with tap water conductivity of about 200uS, the output water limit threshold of 15uS, as recommended by the standards, is reached after about 1500 liters. With a lower quality of the tap water the filters will last for a shorter period and will require more frequent maintenance.

No special maintenance operations are required from the user; the only component that can be user-inspected is the input pre-filter that could get clogged in case of maintenance interventions on the main waterline.

4.2 CLEANING OF THE INPUT PRE-FILTER

- Close off the water supply at the tap.
- Disconnect the power supply module from the mains.
- Switch off the sterilizer and wait for it to cool.
- Remove the elbow fitting of the draining pipe from the PURITY back panel connector.
- Using the special key provided with the sterilizer, remove the filter from the back panel connector (the filter is similar to the ones used on the sterilizer).
- Wash the filter under running water.
- Replace the filter in the connector.
- Connect the draining pipe.
- Open the tap and check for no leaks.
- Connect the power supply module and switch on PURITY.

4.3 DIAGNOSTICS

On a cold power-on of the sterilizer and with the main tank full (Led MAX of clean water tank On) the displays will show the following message at the end of the initial auto-test:

**H2O GOOD or
H2O HARD.**

The message **H2O HARD** indicates the presence of high conductivity water and, consequently, the need for PURITY maintenance; if this message is displayed in the first days after installation, it might be possible that the filters are not fully active yet. Before asking for maintenance, we suggest that the clean water tank be completely emptied, letting PURITY to fill it up again.

If you hear the pump working but no water comes out of PURITY, check that the water tap is open and if necessary perform the cleaning of the pre-filter as described above.

4.4 FILTER KIT REPLACEMENT

With time, the efficiency of the system is reducing more and more, and in the end this will cause the display of the **H2O HARD** message; even though the alarm threshold is passed the rise in conductivity is slow and the water supplied by PURITY can still be of acceptable quality.

WARNING **If PURITY is kept being used for long periods, this may result in damage to the sterilizer or to the tools under sterilization.**

The maintenance intervention by authorized technicians is therefore necessary and cannot be postponed; the filter cartridges must be replaced with a new set. The instructions for the filters replacement are attached to the regeneration kit.

During the annual maintenance of the sterilizer is a good practice to have the quality of the water produced by PURITY checked with an external conductivity meter; the replacement of the filters is recommended when the filters have been used for more than 2000 sterilization cycles.

Take note of the date of filter replacement or of the measured water conductivity value on the SERVICE BOOK attached to the sterilizer's operating manual.

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