



| THERMAL DISINFECTION BY MEANS A BUILT-IN STEAM GENERATOR | | | | | | |
|----------------------------------------------------------|----------|-----------|-------------|---------------|------------------|-------|
| "X" | VOLTAGE | FREQUENCY | RATED POWER | MAXIMUM POWER | ABSORBED CURRENT | FUSES |
| | 200V ~ | 50-60 Hz | 4050 W | 3200 W | 18 A | 20 A |
| | 200V 3~ | | | | 10 A | 16 A |
| | 230V ~ | | | | 15 A | |
| | 230V 3~ | | | | 9 A | |
| | 400V 3N~ | | | | 5 A | |

| | | | | |
|------|-----------|------|------|-----------------------------------|
| 900 | 112402004 | 0 | 1 | C-LINE_AF2.686G - LAVAPADELLE |
| 170 | 119902735 | 0 | 1 | APPL.PIANO SUPERIORE INOX AF2.686 |
| RIF. | CODICE | REV. | Q.TÀ | DESCRIZIONE |

| NOTE: | |
|------------------|--------------|
| PROGETTO/PROJECT | WORKTOP ASAP |
| STANZA/ROOM/RAUM | DIRTY ROOM |

| | | | |
|----------------------|---------------------------|----------------------|------------|
| REFERENCE OFFER NR.: | ORD. CONF. REFERENCE NR.: | FIRST EMISSION DATE: | 12/03/2018 |
|----------------------|---------------------------|----------------------|------------|

OBJECT: **SLUICE ROOM AF2.686/DESTO INSTALLATION - DIM.MM. L600x P605x H900**

| | | | | | |
|--|-------------------------------------------------------------------------------------------------------------|------------|-------|-----------|--------------------|
| | <p>GENERAL TOLLERANCES FOR DIMENTIONES AND ELEMENTS WITHOUT SPECIFICS TOLLERANCES: UNI ISO 2768 - m</p> | PAPER SIZE | A3 | OP.: | C.A. |
| | | SCALE | 1:15 | REV. NR.: | 00 DATE 21/02/2013 |
| | | PAGE | 1 / 1 | DWG NR.: | 180075-00 (CODE) |

LEGENDA. (Marcare con "X" gli allacciamenti in uso)

| X | Pos. | Descrizione |
|----|------|------------------------------------------|
| A | | HOT WATER CONNECTION-3/4"G (TUBE 1/2"G) |
| B | | COLD WATER CONNECTION-3/4"G (TUBE 1/2"G) |
| E1 | | DRAIN TRAP CONNECTION ø110 (A WALL) |
| E2 | | DRAIN TRAP CONNECTION ø110 (A FLOOR) |
| G | | ELECTRICAL CONNECTION BOX |