



THERMAL DISINFECTION BY MEANS A BUILT-IN STEAM GENERATOR						
"X"	VOLTAGE	FREQUENCY	RATED POWER	MAXIMUM POWER	ABSORBED CURRENT	FUSES
	200V ~	50-60 Hz	4050 W	3200 W	18 A	20 A
	200V 3~				10 A	
	230V ~				15 A	16 A
	230V 3~				9 A	
	400V 3N~				5 A	

LEGEND. (Mark with "X" the connections in use)

X	Pos.	Description
A		HOT WATER CONNECTION-3/4"G (TUBE 1/2"G)
B		COLD WATER CONNECTION-3/4"G (TUBE 1/2"G)
E1		DRAIN TRAP CONNECTION DN100 (A WALL)
E2		DRAIN TRAP CONNECTION DN100 (A FLOOR)
F		OMNIPOLAR MAGNETO SWITCH
G		ELECTRICAL CONNECTION BOX
H		EQUIPOTENTIAL CLAMP

NOTE:

REFERENCE OFFER NR.:	ORD. CONF. REFERENCE NR.:	FIRST EMISSION DATE:	01/07/2010
OBJECT: AF2-75B.G-860 - DIM.MM. L750xP605xH900			
	<p>GENERAL TOLERANCES FOR DIMENSIONS AND ELEMENTS WITHOUT SPECIFIC TOLERANCES: UNI ISO 2768 - m</p>	PAPER SIZE	A3
		SCALE	1:12
		PAGE	4 / 5
		OP.:	C.A.
		REV. NR.:	0 DATE
		DWG NR.:	AF2-75B.G-860
		(CODE)	