



THERMAL DISINFECTION BY MEANS A BUILT-IN STEAM GENERATOR

"X"	VOLTAGE	FREQUENCY	RATED POWER	MAXIMUM POWER	ABSORB CURRENT	MAIN SWITCH
	200V ~	50-60 Hz	2700 W	2050 W	9A	16 A
	200V 3~				7A	
	230V ~		11A			
	230V 3~		7A			
	400V 3N~		3550 W	2700 W	7A	

LEGEND. (Mark with "X" the connections in use)

X	Pos	Description
	A	HOT WATER CONNECTION-3/4"G (TUBE 1/2"G)
	B	COLD WATER CONNECTION-3/4"G (TUBE 1/2"G)
	E1	DRAIN TRAP CONNECTION DN100 (A WALL)
	E2	DRAIN TRAP CONNECTION DN100 (A FLOOR)
	F	OMNIPOLAR MAGNETO SWITCH
	G	ELECTRICAL CONNECTION BOX
	H	EQUIPOTENTIAL CLAMP

NOTE:

REFERENCE OFFER NR.:	ORD. CONF. REFERENCE NR.:	FIRST EMISSION DATE:	10/07/2012	
OBJECT: AF2 45PEG - COMPACT - DIM.MM.L450xP500xH1500				
<p>GENERAL TOLLERANCES FOR DIMENTIONES AND ELEMENTS WITHOUT SPECIFICS TOLLERANCES: UNI ISO 2768 - m</p>		PAPER SIZE	A3	
	SCALE	1:10	OP.:	C.A.
	PAGE	4 / 6	REV. NR.:	0 DATE
		DWG NR.:	(CODE) C-AF2 45PEG	