



THERMAL DISINFECTION BY MEANS A BUILT-IN STEAM GENERATOR						
"X"	VOLTAGE	FREQUENCY	RATED POWER	MAXIMUM POWER	ABSORB CURRENT	MAIN SWITCH
	200V ~	50-60 Hz	2700 W	2050 W	9A	16 A
	200V 3~				7A	
	230V ~		3550 W	2700 W	11A	
	230V 3~				7A	
	400V 3N~			7A		

LEGEND. (Mark with "X" the connections in use)

X	Pos	Description
A		HOT WATER CONNECTION-3/4"G (TUBE 1/2"G)
B		COLD WATER CONNECTION-3/4"G (TUBE 1/2"G)
E1		DRAIN TRAP CONNECTION ø110 (A WALL)
E2		DRAIN TRAP CONNECTION ø110 (A FLOW)
F		OMNIPOLAR MAGNETO SWITCH
G		ELECTRICAL CONNECTION BOX
H		EQUIPOTENTIAL CLAMP

NOTE:

REFERENCE OFFER NR.: _____ ORD. CONF. REFERENCE NR.: _____ FIRST EMISSION DATE: 03/07/2012

OBJECT: **AF290BG - PROFESSIONAL - DIM.MM. L900xP600xH1000**



GENERAL TOLLERANCES FOR DIMENTIONES AND ELEMENTS WITHOUT SPECIFICS TOLLERANCES:
UNI ISO 2768 - m



PAPER SIZE A3
SCALE 1:15
PAGE 2 / 5

OP.: C.A.
REV. NR.: 0 DATE
DWG NR.: (CODE) P-AF290BG